


#### **EXPLORATORY RESEARCH: CORE POWERTRAIN MATERIALS**

# DEVELOPMENT OF HIGH TEMPERATURE SAMPLE ENVIRONMENT FOR ADVANCED ALLOY CHARACTERIZATION USING HIGH-ENERGY X-RAY TECHNIQUES

Project ID: mat179



PI: ANDREW CHUANG

Team Members: Peter Kenesei, Jonathon Almer, Dileep Singh X-ray Science Division, Argonne National Laboratory 2020 DOE Vehicle Technology Office Annual Merit Review June 1-4, 2020 Arlington, VA

This presentation does not contain any proprietary, confidential, or otherwise restricted information

# **OVERVIEW**

#### **Timeline**

- Project start: March 2019
- Project end: March 2020
- 90% complete

## **Budget**

FY19 = \$100 K (DOE)

## **Project Partners**

 Applied Materials Division, Argonne National Laboratory

#### **Barriers**

- Performance: Flexible heating device for various material system and specimen geometry. Expand temperature limit up to 1400 °C with high heating and cooling rate (>10 °C /sec) to probe various part of TTT diagram.
- Thermal stability: Device should have minimal thermal load to the surroundings. Requires high thermal stability to minimize motion blur and improve image quality



# RELEVANCE

### Develop characterization capability to support alloy development research

#### Motivation:

- A strong and growing interest to access the most advanced synchrotron X-ray techniques for in-situ material characterization.
- To understand material behavior under in-service conditions requires information at elevated temperature

#### Objectives:

- Develop flexible high temperature sample environment that allows µm stability for high resolution,
   zoom-in/out imaging and diffraction techniques
- Provide improved sample environment accessibility and flexible sample environment, avoid duplicate effort among community, enhance scientific productivity, and enable expansion of in situ testing capabilities.


#### Impact

- Facilitate beamline access for advance powertrain materials research
- Provide the powertrain community with access to a wide range of X-ray techniques at the APS.




# **APPROACH**

## Full-field X-ray Imaging at Elevated Temperature



- Locate region of interest in the bulk sample with regular μtomography (resolution 1~2 μm)
- Use x-ray lens view internal structure of polycrystalline materials such as voids, cracks and inclusions. (resolution ~0.5 μm [now], ~0.1 μm [future])
- Proposed project adds high temperature environment to the high energy x-ray beamline to study material under in-service condition



# **APPROACH**

## Rapid heating/cooling with high thermal stability

Induction based heating

Flexible heating geometry with custom coil

High heating and cooling rate ( > 20°C/sec)

Heat only the specimen, minimize heat load to the environment. Enables µm level imaging.

Use susceptor to increase temperature uniformity


Closed-loop temperature control.


Tempera monitor by thermocouple and thermo

camera

Drive shaft material











## **Technical Accomplishments**

# **HIGH TEMPERATURE 3D IMAGING SETUP**




- Susceptor in quartz tube.
- Whole structure moves down to cover sample.
- Sample is placed at center of susceptor.



Device heating/cooling rate >25 °C/sec

Sample on ceramic(quartz) rod. Glue with high temperature paste







# **Technical Accomplishments**


# TOMOGRAPHY AT HIGH TEMPERATURE

Flush with Ar during experiment





Tomography Image at 900 °C





# PROJECT MILESTONE AND COLLABORATIONS

#### Milestone

Milestone	Description	Status
System design, development and test	<ul> <li>Design and construct the system</li> <li>Test in-house without X-ray</li> </ul>	<ul> <li>Close loop temperature control system finished and tested</li> <li>Reach temperature target of 1350 °C with stability ± 5°C</li> </ul>
Commission the system to work with existing beamline equipment	<ul> <li>Test the system at the beamline</li> <li>Conduct x-ray imaging and diffraction experiment and elevated temperature</li> </ul>	<ul> <li>Conducted two user experiment</li> <li>Collect tomography image at 900°C and diffraction at 1150°C</li> </ul>

#### **Collaborations**

- MINES Saint-Etienne, France (Andras Borbely) In-situ study of recovery and recrystallization of an AlScZr alloy
- Caltech, CA, USA (Benjamin Herren, Katherine Faber) Damage Evolution in Ceramic-Matrix Composite/Environmental Barrier Coating Systems



# PROPOSED FUTURE RESEARCH\*

(\*Any proposed future work is subject to change based on funding levels)

 Add and integrate a compact load frame to the system to characterize engineering material under load at elevated temperature.

# PROJECT SUMMARY

- Rapid heating and cooling cell for high-resolution, high-energy x-ray imaging and diffraction techniques has been developed.
- The system utilizes the concept of induction heating to enable rapid heating and cooling (> 25°C/sec) capability while minimizing the thermal load to the surrounding to achieve ~µm image resolution.
- Proof of concept experiment has been done in late 2019. In a user's in-situ high temperature X-ray experiment, samples (<1g) was heated up to 1250 °C for diffraction experiment and 900 °C for tomography experiment.

